

20 September 2019

## **Medtech Evolution MT2PenCAT QIIP Update Notification**

Medtech Global is pleased to make available for download, the Medtech Evolution MT2PenCAT QIIP release update. This tool helps in mapping the Data elements of Clinical Management system and data elements used in the PENCAT and allow to extract the data in the XML format. This XML file shall be uploaded in the specified path. After which the XML files will be imported in PENCAT tool to generate the required PIP QI Measures.

### **The following are the QI PIP Measures:**

1. Proportion of patient with Smoking status
2. Proportion of patients with Alcohol consumption
3. Proportion of patients with diabetes with no HbA1c result
4. Proportion of patients with diabetes with a current HbA1c
5. Proportion of patients with a weight classification
6. Proportion of patients aged 65 and over who were immunized against Influenza
7. Proportion of patients with diabetes who were immunized against influenza
8. Proportion of patients with COPD who were immunized against influenza
9. Proportion of patients with the necessary risk factors assessed to enable CVD assessment
10. Proportion of female patients with an up-to-date cervical screening
11. Proportion of patients with diabetes with a blood pressure result

## **Important Information – Download and Installation**

**You must have Version 10.3.0 Build 5418 (or above) installed to run the MT2PENCAT update.**

This update MUST be run on the Medtech Evolution Server. All users must be logged out while this update is being run. Please ensure that you have downloaded and read the relevant Release Notes and Installation Instructions provided for this update. Medtech recommend that a suitably qualified IT engineer carries out this upgrade, as a complete Backup is required before the installation of MT2PENCAT update.

**Please refer to the Installation Instructions and Release Notes for more detailed information.**

1. Please visit our website at: <http://www.medtechglobal.com/au/>
2. Select Region: Australia from the Region dropdown list on the top right corner of the screen,
3. Select from the top menu, Support ► Downloads ► Evolution Downloads
4. Locate the update called Version MT2PENCAT Update (September 2019)
5. Review the relevant Release Notes and Installation Instructions
6. Download and run the Update

If you require more information or assistance, please contact Medtech Support:

- Via the Medtech Evolution application [Help ► Contact Support]
- Email: [supportau@medtechglobal.com](mailto:supportau@medtechglobal.com)
- Phone: 1800 148 165 Option 1

Kind Regards,

Customer Care Team